


Rabbit Anti-Human Cd146 Monoclonal Antibody For Fc - P43121

Item Number: CM0008643



Introduction

Rabbit anti-Human CD146 monoclonal antibody for flow cytometry (FC). Validated for human CD146/MUC18. Suitable for studying cell adhesion, melanoma metastasis, and endothelial biology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD146 mAb
Remarks/Alias	MCAM; CD146; MUC18; cell surface glycoprotein MUC18
Species	Human
GeneID (Human)	4162
GeneID	4162
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 24-559 of human CD146 (NP_006491.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P43121
Protein Weight	57kDa/71kDa
Shipping	Ice bag