

Rabbit Anti-Human Cd146 Monoclonal Antibody For Flow Cytometry - P43121

Item Number: CM0012682



Introduction

-conjugated rabbit monoclonal antibody against human CD146 (MUC18/MCAM) for flow cytometry. Validated for human samples; targets a cell adhesion molecule involved in melanoma metastasis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD146 mAb
Remarks/Alias	CD146; MUC18; HEMCAM; METCAM; MeICAM
Species	Human
Gene ID (Human)	4162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-559 of human CD146 (NP_006491.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P43121
Protein Weight	57kDa/72kDa
Shipping	Ice bag