

Rabbit Anti-Human Cd150/Slam Monoclonal Antibody For Fc - Q13291

Item Number: CM0011948



Introduction

-conjugated rabbit monoclonal antibody against human CD150/SLAM, validated for flow cytometry. Targets the extracellular domain of the signaling lymphocytic activation molecule; ideal for immunological studies of T cells, B cells, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD150/SLAM mAb
Remarks/Alias	SLAM; CD150; CDw150
Species	Human
Gene ID (Human)	6504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-237 of human CD150/SLAM (NP_003028.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13291
Protein Weight	37kDa
Shipping	Ice bag