

Rabbit Anti-Human Cd150/Slam Monoclonal Antibody For Fc - Q13291

Item Number: CM0008995



Introduction

-conjugated rabbit monoclonal antibody targeting human CD150/SLAM. Validated for flow cytometry applications. Suitable for immune cell profiling and SLAMF1 signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD150/SLAM mAb
Remarks/Alias	SLAM; CD150; CDw150
Species	Human
GeneID (Human)	6504
GeneID	6504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-237 of human CD150 (NP_003028.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13291
Protein Weight	37kDa
Shipping	Ice bag