

Rabbit Anti-Human Cd150/Slam Monoclonal Antibody For Flow Cytometry - Q13291

Item Number: CM0010562



Introduction

Rabbit monoclonal antibody conjugated to , specific for human CD150/SLAM (SLAMF1), validated for flow cytometry. Targets the signaling lymphocytic activation molecule involved in immune cell regulation and host-pathogen interaction.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD150/SLAM mAb
Remarks/Alias	SLAM; CD150; CDw150
Species	Human
Gene ID (Human)	6504
Gene ID	6504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-237 of human CD150 (NP_003028.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13291
Protein Weight	37kDa
Shipping	Ice bag