

Rabbit Anti-Human Cd152/Ctla-4 Monoclonal Antibody For Flow Cytometry (Intracellular) - P16410

Item Number: CM0005781



Introduction

-conjugated rabbit monoclonal antibody against human CTLA-4 (CD152). For intracellular flow cytometry. Validated in human samples. An immune checkpoint inhibitor target, crucial for T-cell regulation research.

[Learn More](#)

Parameter	Detail
Product Name	Rabbit anti-Human CD152/CTLA-4 mAb
Remarks/Alias	CD; GSE; GRD4; ALPS5; CD152; CTLA-4; IDDM12; CELIAC3
Species	Human
GeneID	1493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-162 of human CD152/CTLA-4 (NP_005205.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P16410
Protein Weight	25kDa
Shipping	Ice bag