


Rabbit Anti-Human Cd157/Bst-1 Mab For Fc - Q10588

Item Number: CM0009617



Introduction

Rabbit anti-Human CD157/BST-1 monoclonal antibody validated for flow cytometry. Targets human CD157/BST-1 (SWISS Q10588), a GPI-anchored ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase involved in calcium signaling and pre-B-cell growth.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD157/BST-1 mAb
Remarks/Alias	CD157; cADPR2; BST1
Species	Human
Gene ID (Human)	683
Gene ID	683
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-292 of human CD157/BST-1 (NP_004325.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q10588
Protein Weight	36kDa/29kDa
Shipping	Ice bag