


## Rabbit Anti-Human Cd157/Bst-1 Mab For Flow Cytometry - Q10588

Item Number: CM0010958



### Introduction

-conjugated rabbit monoclonal antibody targeting human CD157/BST-1. Validated for flow cytometry, reacts with human. Ideal for studying CD157/BST-1 in calcium signaling and pre-B-cell growth.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD157/BST-1 mAb                                                                                       |
| Remarks/Alias    | CD157; cADPR2; BST1                                                                                                     |
| Species          | Human                                                                                                                   |
| GeneID (human)   | 683                                                                                                                     |
| GeneID           | 683                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 29-292 of human CD157/BST-1 (NP_004325.2) |
| Source           | Rabbit                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                   |
| Application      | FC                                                                                                                      |
| Cross Reactivity | Human                                                                                                                   |
| SWISS            | Q10588                                                                                                                  |
| Protein Weight   | 36kDa/29kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                                 |