

## Rabbit Anti-Human Cd158 (Kir2DL1/S1/S2/S3/S5) Monoclonal Antibody For Flow Cytometry - Q14954 / Q14953 / Q14952 / P43631 / P43626

Item Number: CM0011072

### Introduction

High-quality -conjugated rabbit monoclonal antibody targeting human CD158 (KIR2DL1/S1/S2/S3/S5) for flow cytometry. Validated for human samples, this antibody enables precise NK cell receptor analysis in immunology research.

[Learn More](#)

| Parameter        | Value                                                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD158 (KIR2DL1/S1/S2/S3/S5) mAb                                                                                        |
| Remarks/Alias    | NKAT9; CD158G CD158A NKAT1; CD158H; CD158J; NKAT5                                                                                        |
| Species          | Human                                                                                                                                    |
| GeneID (Human)   | 3810/3808/3806/3802/100132285                                                                                                            |
| GeneID           | 3810/3808/3806/3802/100132285                                                                                                            |
| Immunogen        | Recombinant protein   Recombinant fusion protein containing a sequence corresponding to amino acids 22-245 of human CD158 (NP_055327.1). |
| Source           | Rabbit                                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                                    |
| Application      | FC                                                                                                                                       |
| Cross Reactivity | Human                                                                                                                                    |
| SWISS            | Q14954/Q14953/Q14952/P43631/P43626                                                                                                       |
| Protein Weight   | 34kDa                                                                                                                                    |
| Shipping         | Ice bag                                                                                                                                  |