

Rabbit Anti-Human Cd158A/Kir2DL1 Monoclonal Antibody For Flow Cytometry - P43626

Item Number: CM0004926



Introduction

-conjugated rabbit monoclonal antibody targeting human CD158a/KIR2DL1, validated for flow cytometry. Ideal for NK cell receptor studies and immune profiling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD158a/KIR2DL1 mAb
Remarks/Alias	NKAT; NKAT1; p58.1; CD158A; KIR221; NKAT-1; KIR-K64; KIR2DL3
Species	Human
GeneID (Human)	3802
GeneID	3802
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-245 of human CD158a/KIR2DL1 (NP_055033.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P43626
Protein Weight	38kDa/41kDa
Shipping	Ice bag