

Rabbit Anti-Human Cd161 Monoclonal Antibody For Flow Cytometry - Q12918

Item Number: CM0004114



Introduction

-conjugated rabbit monoclonal antibody against human CD161 (KLRB1, NKR-P1A) for flow cytometry. Specific for human samples, detects a 25 kDa membrane receptor involved in NK cell inhibition.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD161 mAb
Remarks/Alias	NKR; CD161; CLEC5B; NKR-P1; NKR-P1A; hNKR-P1A; KLRB1
Species	Human
Gene ID (Human)	3820
Gene ID	3820
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-225 of human CD161 (NP_002249.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q12918
Protein Weight	25kDa
Shipping	Ice bag