

Rabbit Anti-Human Cd162/Psgl-1 Monoclonal Antibody For Flow Cytometry - Q14242

Item Number: CM0003998



Introduction

-conjugated rabbit monoclonal antibody targeting human CD162/PSGL-1, a selectin ligand. Validated for flow cytometry (FC) on human leukocytes. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD162/PSGL-1 mAb
Remarks/Alias	CLA; CD162; PSGL1; PSGL-1
Species	Human
GeneID (Human)	6404
GeneID	6404
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-295 of human CD162/PSGL-1 (NP_002997.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q14242
Protein Weight	43kDa/44kDa
Shipping	Ice bag