


Rabbit Anti-Human Cd166/Alcam Monoclonal Antibody For Fc - Q13740

Item Number: CM0008976



Introduction

High-quality -conjugated rabbit monoclonal antibody against human CD166/ALCAM. Validated for flow cytometry (FC). Ideal for T-cell activation, cell adhesion, and hematopoietic stem cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD166/ALCAM mAb
Remarks/Alias	ALCAM; CD166; MEMD; CD166 antigen
Species	Human
Gene ID (Human)	214
Gene ID	214
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13740
Protein Weight	15kDa/33kDa/63kDa/65kDa
Shipping	Ice bag