


Rabbit Anti-Human Cd166/Alcam Monoclonal Antibody For Flow Cytometry - Q13740

Item Number: CM0009334



Introduction

-conjugated rabbit monoclonal antibody specific for human CD166/ALCAM (UniProt Q13740), designed for flow cytometry. CD166 mediates cell adhesion, T-cell activation, and hematopoietic stem cell engraftment. Ideal for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD166/ALCAM mAb
Remarks/Alias	MEMD; CD166
Species	Human
GeneID (Human)	214
GeneID	214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-526 of human CD166/ALCAM (NP_001618.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13740
Protein Weight	15kDa/33kDa/63kDa/65kDa
Shipping	Ice bag