

## Rabbit Anti-Human Cd16B Mab For Fc - 075015

Item Number: CM0008051



### Introduction

-conjugated rabbit monoclonal antibody against human CD16b (FCGR3B). Designed for flow cytometry, detects low-affinity IgG Fc receptor on neutrophils. Useful for immunology and immune complex research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD16b mAb                                                                                       |
| Remarks/Alias    | CD16, CD16-I, CD16b, FCG3, FCGR3, FCGR3B, FCR-10, FCRIII, FCRIIIb                                                 |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 2215                                                                                                              |
| GeneID           | 2215                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 20-208 of human CD16b (NP_000561.3) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Human                                                                                                             |
| SWISS            | 075015                                                                                                            |
| Protein Weight   | 26KDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |