

Rabbit Anti-Human Cd16B Monoclonal Antibody For Flow Cytometry - O75015

Item Number: CM0009097



Introduction

-conjugated rabbit monoclonal antibody targeting human CD16b (FcγRIIIb) for flow cytometry. Low affinity IgG Fc receptor expressed on neutrophils; suitable for immune complex studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD16b mAb
Remarks/Alias	CD16, CD16-I, CD16b, FCG3, FCGR3, FCGR3B, FCR-10, FCRIII, FCRIIIb
Species	Human
Gene ID	2215
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-208 of human CD16b (NP_000561.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75015
Protein Weight	26KDa
Shipping	Ice bag