


Rabbit Anti-Human Cd170/Siglec-5 Mab For Fc - O15389

Item Number: CM0011369



Introduction

-conjugated rabbit monoclonal antibody targeting human CD170/Siglec-5, validated for flow cytometry. Detects sialic acid-binding Ig-like lectin 5, an adhesion molecule on myeloid cells. Suitable for immunophenotyping studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD170/Siglec-5 mAb
Remarks/Alias	CD170, OBBP2, CD33L2, OB-BP2, SIGLEC-5
Species	Human
GeneID (Human)	8778
GeneID	8778
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 17-434 of human CD170/Siglec-5 (NP_003821.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O15389
Protein Weight	61kDa
Shipping	Ice bag