

Rabbit Anti-Human Cd171/L1Cam Monoclonal Antibody For Flow Cytometry - P32004

Item Number: CM0008875



Introduction

-conjugated rabbit monoclonal antibody targeting human CD171/L1CAM (neural cell adhesion molecule L1). Validated for flow cytometry, suitable for neuroscience, neuronal migration, and synaptic plasticity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD171/L1CAM mAb
Remarks/Alias	S10; HSAS; HYCX; MASA; MIC5; SPG1; CAML1; CD171; HSAS1; N-CAML1; NCAM-L1; N-CAM-L1
Species	Human
Gene ID (Human)	3897
Gene ID	3897
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-1120 of human CD171/L1CAM (NP_000416.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32004
Protein Weight	140kDa/139kDa
Shipping	Ice bag