

Abflo 647 Rabbit Anti-Human Cd172A/Sirpα Monoclonal Antibody For Flow Cytometry - P78324

Item Number: CM0010430



Introduction

ABflo 647-conjugated rabbit monoclonal antibody against human CD172a/SIRPα (UniProt P78324). Validated for flow cytometry. Ideal for studying CD47-SIRPα signaling, phagocytosis regulation, and immune checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD172a/SIRPα mAb
Remarks/Alias	BIT; MFR; P84; SIRP; MYD-1; SHPS1; CD172A; PTPNS1
Species	Human
Gene ID	140885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-370 of human CD172a/SIRPα (NP_542970.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P78324
Protein Weight	55kDa
Shipping	Ice bag