

Rabbit Anti-Human Cd172B Monoclonal Antibody For Flow Cytometry - O00241

Item Number: CM0009457



Introduction

-conjugated rabbit monoclonal antibody targeting human CD172b (SIRP-beta-1). Validated for flow cytometry, reacts with human. Recombinant immunogen, molecular weight 43 kDa. Ideal for studying myeloid cell signaling and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD172b mAb
Remarks/Alias	CD172b; SIRP-BETA-1
Species	Human
GeneID	10326
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-371 of human CD172b (NP_006056.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O00241
Protein Weight	43kDa
Shipping	Ice bag