

Rabbit Anti-Human Cd177 Monoclonal Antibody () For Flow Cytometry - Q8N6Q3

Item Number: CM0010325

Introduction

Rabbit monoclonal antibody against human CD177, conjugated to for flow cytometry. Recognizes a neutrophil-specific GPI-anchored receptor involved in activation and migration. Suitable for immunology, inflammation, and hematological research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD177 mAb
Remarks/Alias	CD177; HNA-2a; HNA2A; NB1; NB1 GP; PRV-1; PRV1; NB1GP; CD177 antigen
Species	Human
GeneID (Human)	57126
GeneID	57126
Immunogen	Recombinant fusion protein containing amino acids 22-407 of human CD177 (NP_065139.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human
SWISS	Q8N6Q3
Protein Weight	15kDa/27kDa/46kDa
Shipping	Ice bag