

Rabbit Anti-Human Cd185/Cxcr5 Monoclonal Antibody For Fc - P32302

Item Number: CM0005585

Introduction

-conjugated rabbit monoclonal antibody against human CD185/CXCR5. Designed for flow cytometry, useful for B-cell migration and Burkitt lymphoma research. Confirmed cross-reactivity with human.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD185/CXCR5 mAb
Remarks/Alias	BLR1; CD185; MDR15
Species	Human
GenelD (human)	643
GenelD	643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-47aa of Human CD185/CXCR5 (NP_033981.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32302
Protein Weight	42kDa
Shipping	Ice bag