

Rabbit Anti-Human Cd185/Cxcr5 Monoclonal Antibody For Flow Cytometry - P32302

Item Number: CM0009055



Introduction

Rabbit anti-Human CD185/CXCR5 monoclonal antibody validated for flow cytometry, ideal for B-cell migration and Burkitt lymphoma research. One-stop IVD raw material sourcing from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD185/CXCR5 mAb
Remarks/Alias	BLR1; CD185; MDR15
Species	Human
GeneID (human)	643
GeneID	643
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-47 of Human CD185/CXCR5 (NP_001707.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32302
Protein Weight	42kDa
Shipping	Ice bag