


Rabbit Anti-Human Cd197/Ccr7 Monoclonal Antibody For Flow Cytometry - P32248

Item Number: CM0007918



Introduction

-conjugated rabbit monoclonal antibody targeting human CD197/CCR7, validated for flow cytometry. Detects a 43kDa GPCR involved in lymphocyte homing, ideal for immunology and EBV research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD197/CCR7 mAb
Remarks/Alias	BLR2; EB11; CCR-7; CD197; CDw197; CMKBR7; CC-CCR-7
Species	Human
Gene ID	1236
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-59 of human CD197/CCR7 (NP_001829.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32248
Protein Weight	43kDa
Shipping	Ice bag