


Rabbit Anti-Human Cd1B Monoclonal Antibody For Fc - P29016

Item Number: CM0011581



Introduction

-conjugated rabbit monoclonal antibody against human CD1b (SWISS P29016) for flow cytometry. Ideal for studying lipid antigen presentation and T-cell responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD1B mAb
Remarks/Alias	R1; CD1; CD1A
Species	Human
Gene ID (Human)	910
Gene ID	910
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-303 of human CD1B (NP_001755.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P29016
Protein Weight	37kDa/30kDa
Shipping	Ice bag