

Rabbit Anti-Human Cd1C Monoclonal Antibody For Fc - P29017

Item Number: CM0012411



Introduction

-conjugated rabbit anti-human CD1c monoclonal antibody, validated for flow cytometry. Targets the lipid antigen-presenting molecule CD1c (UniProt P29017). Ideal for studying T-cell immunology, antigen presentation, and lymphoid malignancies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD1c mAb
Remarks/Alias	R7; CD1; CD1A; BDCA1
Species	Human
GeneID (human)	911
GeneID	911
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-298 of human CD1c (NP_001756.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P29017
Protein Weight	38kDa
Shipping	Ice bag