


Rabbit Anti-Human Cd1D Mab - P15813

Item Number: CM0009531



Introduction

-conjugated rabbit monoclonal antibody targeting human CD1d, a key antigen-presenting glycoprotein for glycolipids. Validated for Flow Cytometry. Ideal for NKT cell biology and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD1d mAb
Remarks/Alias	R3; CD1A; R3G1
Species	Human
GeneID (Human)	912
GeneID	912
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-301 of human CD1d (NP_001757.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15813
Protein Weight	38kDa
Shipping	Ice bag