

Rabbit Anti-Human Cd1D Monoclonal Antibody For Flow Cytometry - P15813

Item Number: CM0011989



Introduction

-conjugated rabbit monoclonal antibody targeting human CD1d (P15813). Validated for flow cytometry to study glycolipid presentation and NKT cell biology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD1d mAb
Remarks/Alias	R3; CD1A; R3G1
Species	Human
Gene ID (Human)	912
Gene ID	912
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to aminoacids 20-301 of human CD1d (NP_001757.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15813
Protein Weight	38kDa
Shipping	Ice bag