

Rabbit Anti-Human Cd1D Monoclonal Antibody For Flow Cytometry - P15813

Item Number: CM0005991

Introduction

-conjugated rabbit monoclonal antibody against human CD1d (clone R3G1). Validated for flow cytometry, targeting the glycolipid-presenting receptor (UniProt P15813). Ideal for NKT cell and innate immunity studies.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD1d mAb
Remarks/Alias	R3, CD1A, R3G1
Species	Human
GeneID (Human)	912
GeneID	912
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-301 of human CD1d (NP_001757.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15813
Protein Weight	38kDa
Shipping	Ice bag