


Rabbit Anti-Human Cd2 Mab - P06729

Item Number: CM0012209



Introduction

-conjugated rabbit monoclonal antibody targeting human CD2 (LFA-2, SRBC), validated for flow cytometry. Ideal for T-cell and NK cell immunophenotyping and adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD2 mAb
Remarks/Alias	T11; SRBC; LFA-2
Species	Human
Gene ID (Human)	914
Gene ID	914
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag