


Rabbit Anti-Human Cd2 Monoclonal Antibody For Flow Cytometry - P06729

Item Number: CM0010833



Introduction

-labeled rabbit monoclonal antibody targeting human CD2 (T11, LFA-2), validated for flow cytometry. Recognizes a 39 kDa T-cell surface glycoprotein involved in adhesion and signaling. Ideal for lymphocyte subset analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD2 mAb
Remarks/Alias	T11; SRBC; LFA-2
Species	Human
Gene ID	914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-209 of human CD2 (NP_001758.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag