

Rabbit Anti-Human Cd2 Monoclonal Antibody For Flow Cytometry - P06729

Item Number: CM0008269



Introduction

-conjugated rabbit monoclonal antibody targeting human CD2, validated for flow cytometry detection. Suitable for T-cell and NK cell research applications. IVD raw material for diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD2 mAb
Remarks/Alias	CD2; LFA-2; SRBC; T11; CD2 molecule
Species	Human
Gene ID	914 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-209 of human CD2 (NP_001758.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag