

## Rabbit Anti-Human Cd200/Ox2 Monoclonal Antibody For Flow Cytometry - P41217

Item Number: CM0008937

### Introduction

-conjugated rabbit monoclonal antibody targeting human CD200/OX2. Validated for flow cytometry applications on human samples. Useful for investigating T-cell costimulation and myeloid cell activity. Research-use only.

[Learn More](#)

| Parameter        | Value                                        |
|------------------|----------------------------------------------|
| Product Name     | Rabbit anti-Human CD200/OX2 mAb              |
| Remarks/Alias    | CD200; MOX1; MOX2; MRC; OX-2; CD200 molecule |
| Species          | Human                                        |
| Gene ID          | 4345                                         |
| Immunogen        | Recombinant protein                          |
| Source           | Rabbit                                       |
| Category         | Monoclonal Antibodies                        |
| Application      | FC                                           |
| Cross Reactivity | Human                                        |
| SWISS            | P41217                                       |
| Protein Weight   | 30kDa/31kDa/32kDa                            |
| Shipping         | Ice bag                                      |