


Rabbit Anti-Human Cd200R Monoclonal Antibody For Flow Cytometry - Q8Td46

Item Number: CM0009261



Introduction

-conjugated rabbit monoclonal antibody targeting human CD200R (OX2R). Designed for flow cytometry (FC) applications. Specific for human CD200 receptor 1, an inhibitory immune checkpoint receptor involved in regulating inflammation. Ideal for immunophenotyping and immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD200R mAb
Remarks/Alias	OX2R; MOX2R; CD200R; HCRTR2
Species	Human
Gene ID (Human)	131450
Gene ID	131450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-240 of human CD200R (NP_740750.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TD46
Protein Weight	19kDa/21kDa/37kDa/39kDa
Shipping	Ice bag