


Rabbit Anti-Human Cd200R Monoclonal Antibody For Flow Cytometry - Q8Td46

Item Number: CM0009025



Introduction

Rabbit monoclonal antibody against human CD200R, conjugated to , validated for flow cytometry. Detects CD200 receptor expression on immune cells in immuno-oncology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD200R mAb
Remarks/Alias	OX2R; MOX2R; CD200R; HCRTR2
Species	Human
Gene ID (Human)	131450
Gene ID	131450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-240 of human CD200R (NP_740750.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TD46
Protein Weight	19kDa/21kDa/37kDa/39kDa
Shipping	Ice bag