


Rabbit Anti-Human Cd201/Epcr Monoclonal Antibody For Fc - Q9Unn8

Item Number: CM0011356



Introduction

-conjugated rabbit monoclonal antibody targeting human CD201/EPCR (EPCR). Validated for flow cytometry, it binds the endothelial protein C receptor involved in coagulation. Ideal for IVD assay development and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD201/EPCR mAb
Remarks/Alias	PROCR; CCCA; CCD41; EPCR; protein C receptor
Species	Human
Gene ID	10544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-210 of human CD201/EPCR (NP_006395.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UNN8
Protein Weight	26kDa
Shipping	Ice bag