


## Rabbit Anti-Human Cd201/Epcr Monoclonal Antibody For Flow Cytometry - Q9Unn8

Item Number: CM0008301



### Introduction

-conjugated rabbit monoclonal antibody targeting human CD201/EPCR, validated for flow cytometry. Ideal for studying the protein C anticoagulant pathway and endothelial cell biology in hemostasis research.

[Learn More](#)

| Parameter        | Value                                                                                                                                       |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD201/EPCR mAb                                                                                                            |
| Remarks/Alias    | PROCR; CCCA; CCD41; EPCR; protein C receptor                                                                                                |
| Species          | Human                                                                                                                                       |
| GeneID (Human)   | 10544                                                                                                                                       |
| GeneID           | 10544                                                                                                                                       |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-210 of human CD201/EPCR (NP_006395.2). |
| Source           | Rabbit                                                                                                                                      |
| Category         | Monoclonal Antibodies                                                                                                                       |
| Application      | FC                                                                                                                                          |
| Cross Reactivity | Human                                                                                                                                       |
| SWISS            | Q9UNN8                                                                                                                                      |
| Protein Weight   | 26kDa                                                                                                                                       |
| Shipping         | Ice bag                                                                                                                                     |