


Rabbit Anti-Human Cd201/Epcr Monoclonal Antibody For Flow Cytometry - Q9Unn8

Item Number: CM0004546



Introduction

-conjugated rabbit monoclonal antibody targeting human CD201/EPCR, validated for flow cytometry. Ideal for studying protein C anticoagulant pathway in endothelial cells, coagulation research, and hemostasis.

[Learn More](#)

Parameter	Description
Product Name	Rabbit anti-Human CD201/EPCR mAb
Remarks/Alias	PROCR; CCCA; CCD41; EPCR; protein C receptor
Species	Human
GeneID (Human)	10544
GeneID	10544
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 18-210 of human CD201/EPCR (NP_006395.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UNN8
Protein Weight	27kDa
Shipping	Ice bag