


Rabbit Anti-Human Cd204 Mab For Fc - P21757

Item Number: CM0009287



Introduction

-conjugated Rabbit anti-Human CD204 monoclonal antibody validated for flow cytometry (FC). Targets human macrophage scavenger receptor MSR1 (CD204), a key receptor in modified LDL uptake and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD204 mAb
Remarks/Alias	SRA; SR-A; CD204; SR-AI; phSR1; phSR2; SCARA1; SR-AII; SR-AIII
Species	Human
Gene ID (Human)	4481
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 77-451 of human CD204 (NP_619729.1).
Source	Rabbit
Category	Monoclonal Antibody
Application	FC
Cross Reactivity	Human
UniProt Accession	P21757
Protein Weight	40kDa/43kDa/50kDa
Shipping	Ice bag