


Rabbit Anti-Human Cd204 Monoclonal Antibody For Flow Cytometry - P21757

Item Number: CM0010382



Introduction

Conjugated monoclonal antibody against human CD204 (MSR1/SCARA1) for flow cytometry. Validated in human samples, recognizing the scavenger receptor involved in modified LDL uptake and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD204 mAb
Remarks/Alias	SRA; SR-A; CD204; SR-AI; phSR1; phSR2; SCARA1; SR-AII; SR-AIII
Species	Human
GeneID (Human)	4481
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-451 of human CD204 (NP_619729.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21757
Protein Weight	40kDa/43kDa/50kDa
Shipping	Ice bag