


Rabbit Anti-Human Cd206/Mmr Monoclonal Antibody For Fc - P22897

Item Number: CM0011747



Introduction

-conjugated rabbit monoclonal antibody targeting human CD206/MMR for flow cytometry. Designed against recombinant extracellular domain, suitable for macrophage mannose receptor detection in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD206/MMR mAb
Remarks/Alias	MMR; hMR; CD206; MRC1L1; CLEC13D; CLEC13DL; bA541i19.1
Species	Human
GeneID (Human)	4360
GeneID	4360
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-1389 of human Mannose Receptor/CD206 (NP_002429.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22897
Protein Weight	166kDa
Shipping	Ice bag