


Rabbit Anti-Human Cd206/Mmr Monoclonal Antibody For Fc - P22897

Item Number: CM0008122



Introduction

-conjugated Rabbit anti-Human CD206/MMR monoclonal antibody for flow cytometry. Targets human macrophage mannose receptor 1, validated for FC. Suitable for immuno-oncology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD206/MMR mAb
Remarks/Alias	MMR; hMR; CD206; MRC1L1; CLEC13D; CLEC13DL; bA541I19.1
Species	Human
Gene ID	4360
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-1389 of human Mannose Receptor/CD206 (NP_002429.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22897
Protein Weight	166kDa
Shipping	Ice bag