


Rabbit Anti-Human Cd213A1/IL-13R α 1 Monoclonal Antibody For Fc - P78552

Item Number: CM0008290



Introduction

-conjugated rabbit monoclonal antibody targeting human CD213a1/IL-13R α 1, validated for flow cytometry. Binds the extracellular domain, suitable for immunology and cytokine signaling research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD213a1/IL-13R α 1 mAb
Remarks/Alias	IL13RA1; CD213A1; CT19; IL-13Ra; NR4; interleukin-13 receptor subunit alpha-1
Species	Human
Gene ID (Human)	3597
Gene ID	3597
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-343 of human CD213a1/IL-13R α 1 (NP_001551.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P78552
Protein Weight	31kDa/48kDa
Shipping	Ice bag