

Rabbit Anti-Human Cd215/IL-15R Alpha Monoclonal Antibody For Fc - Q13261

Item Number: CM0003997



Introduction

High-affinity -conjugated rabbit monoclonal antibody targeting human CD215/IL-15 receptor alpha for flow cytometry. Validated for human samples, ideal for immune cell receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD215/IL-15R alpha mAb
Remarks/Alias	CD215
Species	Human
Gene ID (Human)	3601
Gene ID	3601
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-205 of human CD215/IL-15R alpha (NP_002180.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13261
Protein Weight	28kDa
Shipping	Ice bag