

Rabbit Anti-Human Cd215/IL-15R Alpha Monoclonal Antibody For Flow Cytometry - Q13261

Item Number: CM0010342



Introduction

-conjugated rabbit monoclonal antibody against human CD215/IL-15R alpha (IL15RA). Validated for flow cytometry. Recognizes human receptor with a protein weight of 15-28 kDa, suitable for immune signaling and receptor studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD215/IL-15R alpha mAb
Remarks/Alias	IL15RA; CD215; interleukin-15 receptor subunit alpha
Species	Human
GeneID	3601
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-205 of human CD215/IL-15R alpha (NP_002180.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13261
Protein Weight	15-28kDa
Shipping	Ice bag