

## Rabbit Anti-Human Cd22 Monoclonal Antibody For Flow Cytometry - P20273

Item Number: CM0006086



### Introduction

-conjugated rabbit anti-human CD22 monoclonal antibody for flow cytometry. Reacts with human B-cell receptor CD22 (95 kDa), a member of the Siglec family involved in B-cell signaling and adhesion. Ideal for B-cell research.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD22 mAb                                                                                                            |
| Remarks/Alias    | CD22; SIGLEC-2; SIGLEC2; CD22 molecule                                                                                                |
| Species          | Human                                                                                                                                 |
| GeneID (Human)   | 933                                                                                                                                   |
| GeneID           | 933                                                                                                                                   |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-687 of human CD22 (NP_001762.2). |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | FC                                                                                                                                    |
| Cross Reactivity | Human                                                                                                                                 |
| SWISS            | P20273                                                                                                                                |
| Protein Weight   | 95kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |