


Rabbit Anti-Human Cd22 Monoclonal Antibody For Flow Cytometry - P20273

Item Number: CM0005724



Introduction

-conjugated rabbit monoclonal antibody targeting human CD22 (Siglec-2). Validated for flow cytometry on B-lymphocytes. Recombinant immunogen (aa 20-687). Ideal for immunophenotyping and B-cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD22 mAb
Remarks/Alias	SIGLEC2; SIGLEC-2
Species	Human
Gene ID (Human)	933
Gene ID	933
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-687 of human CD22 (NP_001762.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20273
Protein Weight	95kDa
Shipping	Ice bag