

Rabbit Anti-Human Cd225/Ifitm1 Monoclonal Antibody For Fc (Intra) - P13164

Item Number: CM0009556

Introduction

-conjugated rabbit monoclonal antibody against human CD225/IFITM1 for intracellular flow cytometry. Ideal for antiviral and cancer research; reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD225/IFITM1 mAb
Remarks/Alias	9-27; CD225; IFI17; LEU13; DSPA2a
Species	Human
Gene ID	8519
Immunogen	Recombinant fusion protein containing a sequence within amino acids 1-80 of human CD225/IFITM1 (NP_003632.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P13164
Protein Weight	14kDa
Shipping	Ice bag