


Rabbit Anti-Human Cd226/Dnam-1 Mab For Fc - Q15762

Item Number: CM0012437



Introduction

-conjugated rabbit monoclonal antibody targeting human CD226/DNAM-1, a costimulatory receptor on T and NK cells. Validated for flow cytometry; ideal for studying immune cell adhesion, cytotoxicity, and lymphocyte signaling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD226/DNAM-1 mAb
Remarks/Alias	PTA1; DNAM1; DNAM-1; TLISA1
Species	Human
Gene ID (Human)	10666
Gene ID	10666
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-247 of human CD226/DNAM-1 (NP_006557.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q15762
Protein Weight	39kDa
Shipping	Ice bag