


Rabbit Anti-Human Cd226/Dnam-1 Mab - Q15762

Item Number: CM0007960



Introduction

-conjugated rabbit monoclonal antibody against human CD226/DNAM-1, validated for flow cytometry. Suitable for studying immune cell adhesion, NK and T-cell cytotoxicity, and tumor immunosurveillance.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD226/DNAM-1 mAb
Remarks/Alias	PTA1; DNAM1; DNAM-1; TLIa1
Species	Human
Gene ID (Human)	10666
Gene ID	10666
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-247 of human CD226/DNAM-1 (NP_006557.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q15762
Protein Weight	39kDa
Shipping	Ice bag